

# Flag Nanobody Agarose Beads

## (FNA001)

### Product Overview

Flag Nanobody Agarose Beads are prepared by conjugating anti-Flag-tag (DYKDDDDK) nanobodies to agarose beads, capable of binding Flag-tagged fusion proteins. This product lacks the light and heavy chains of conventional antibodies, with a dissociation constant at the nM-pM level. 20  $\mu$ L of bead suspension can bind 10–15  $\mu$ g of Flag-tagged fusion proteins. It features high specificity, low background, high affinity, high binding capacity, and ease of use, suitable for immunoprecipitation (IP)/co-immunoprecipitation (CoIP), chromatin immunoprecipitation (ChIP)/RNA-binding protein immunoprecipitation (RIP), enzyme activity assays, mass spectrometry analysis, and other experiments.

### Components

| Product Name                | Cat. No.  | Volume      |
|-----------------------------|-----------|-------------|
| Flag Nanobody Agarose Beads | FNA001-01 | 500 $\mu$ L |
|                             | FNA001-02 | 1 mL        |

### Storage Conditions

Stable for 1 year at 4°C or 2 years at -20°C. Shipped with ice packs.

### Precautions

- For your safety, please wear a lab coat and disposable gloves. This product is for **research use only**.
- For long-term storage, keep at -20°C. Avoid centrifugation, drying, and repeated freeze-thaw cycles.
- The storage buffer is PBS containing 20% ethanol and 0.03% NaN<sub>3</sub>.

### Experimental Protocol

#### 1. Cell Collection:

Use approximately 10<sup>6</sup>–10<sup>7</sup> cells expressing Flag-tagged fusion proteins per immunoprecipitation reaction; adjust cell number based on expression level. Aspirate the medium, wash cells twice with pre-chilled 1×PBS, harvest adherent cells using a cell scraper or trypsin digestion, transfer to a centrifuge tube, centrifuge at 500 g for 3 min at 4°C, and discard the supernatant.

#### 2. Cell Lysis:

- 1) Add protease inhibitors to lysis buffer, resuspend cells in 1.0 mL pre-chilled lysis buffer. For membrane proteins or nuclear/chromatin proteins, use RIPA lysis buffer or standard IP lysis buffer combined with sonication.

- 2) Incubate on ice for 30 min, mixing thoroughly every 10 min or rotating on a rotary mixer at 4°C.
- 3) Centrifuge at 20,000 g for 15 min at 4°C, transfer the supernatant (cell lysate) to a new pre-chilled tube, and discard the pellet.

**Note:** The cell lysate can be stored long-term at -80°C.

**3. Bead Equilibration (Optional):** Mix Flag Nanobody Agarose Beads thoroughly, transfer 20 µL to 500 µL pre-chilled lysis buffer, separate beads on a magnetic stand until the supernatant is clear, and discard the supernatant (tip: cut off the end of the pipette tip when aspirating 20 µL bead suspension).

#### 4. Protein Binding:

- 1) Add cell lysate to the equilibrated Flag Nanobody Agarose Beads (or add 20 µL beads directly to the lysate if Step 3 is skipped). Incubate on a rotary mixer at 4°C for 1–3 h. If needed, save 50 µL lysate as input for Western blot analysis.
- 2) Centrifuge at 2500 g for 3 min at 4°C and discard the supernatant.

**5. Bead Washing:** Resuspend beads in 1.0 mL pre-chilled lysis buffer, wash on a rotary mixer at 4°C for 5 min, centrifuge at 2500 g for 3 min at 4°C, discard the supernatant, and repeat washing twice.

#### 6. Protein Elution (Choose method based on downstream application):

- 1) **SDS-PAGE Sample Buffer Elution:** Resuspend beads in 20–50 µL 2×SDS-sample buffer. Denature thoroughly by heating at 95°C for 10 min to release immunoprecipitated complexes from beads, centrifuge at 2500 g for 3 min at 4°C, and collect the supernatant for SDS-PAGE and Western blot analysis.
- 2) **Acidic Elution:** Add 50 µL 0.2 M glycine (pH 2.5) to elute bound proteins, incubate for 30 s with constant mixing, centrifuge, transfer supernatant to a new tube, and immediately add 5 µL 1.0 M Tris (pH 10.4) to neutralize acidic glycine.

**Note:** Repeat this step to improve elution efficiency.

- 3) **Flag Fusion Protein Elution:** Add excess Flag-peptide for competitive elution.