

HS Taq DNA Polymerase

(M011)

Product Description

HS Taq DNA Polymerase is a **chemically modified Taq DNA Polymerase** that is fully inactive at room temperature. Its activity is only released after heating at 95 °C, effectively preventing non-specific amplification and primer-dimer formation during reaction setup and ramp-up.

Compared with antibody-based hot-start Taq enzymes, HS Taq DNA Polymerase provides **more complete inactivation** and higher stringency. Compared with other chemically modified hot-start Taq enzymes, it requires only **5 minutes for activation** and is fully compatible with standard PCR protocols.

Paired with an optimized buffer system, HS Taq DNA Polymerase minimizes non-specific amplification and primer dimers, delivering exceptional sensitivity and specificity. It is ideal for amplifying low-copy genes from complex templates.

二、Components

Product Name	Cat. No.	Volume
HS Taq DNA Polymerase	M011	5 U/μL, 250 U

Definition of Activity Unit

One unit (U) is defined as the amount of enzyme that incorporates **10 nmol of total nucleotides** into acid-insoluble material within **30 minutes at 74 °C** using activated salmon sperm DNA as the template/primer.

Storage Conditions

Storage: **-30 ~ -15 °C**

Shipping: ≤ 0 °C (ice pack)

Avoid repeated freeze-thaw cycles.

Applications

This product is widely used for DNA amplification from animals, plants, microorganisms, and other sources.

Recommended Uses

Featuring 5'→3' exonuclease activity, high specificity, and high sensitivity, this hot-start enzyme is suitable for:

1. Quantitative real-time PCR (qPCR)
2. cDNA library construction

3. Mutant analysis and construction
4. Genetic disease diagnosis
5. Forensic identification

Primer Design Guidelines

1. The 3'-terminal base of primers is preferably **G or C**.
2. Avoid consecutive mismatches in the last 8 bases at the 3' end.
3. Avoid hairpin structures at the 3' end.
4. T_m values of forward and reverse primers should differ by ≤ 1 °C; optimal T_m: **55 ~ 65 °C** (calculated using Primer Premier 5 recommended).
5. Non-complementary overhangs should **not** be included in T_m calculation.
6. GC content: **40% ~ 60%**.
7. Uniform distribution of A, G, C, T; avoid GC- or AT-rich regions.
8. Avoid complementary sequences of ≥ 5 bases within a primer or between primers; avoid ≥ 3 complementary bases at the 3' ends of primer pairs.
9. Verify primer specificity using **NCBI BLAST** to prevent non-specific amplification.

1. Reaction System (Recommended to prepare on ice)

Component	20 μ L Reaction	50 μ L Reaction	Final Concentration
10 $\times$ Taq Buffer	2 μ L	5 μ L	1 $\times$
Primer 1 (10 μ M)	0.4 μ L	1 μ L	0.2 μ M each
Primer 2 (10 μ M)	0.4 μ L	1 μ L	0.2 μ M each
dNTP Mix (10 mM each)	0.4 μ L	1 μ L	0.2 μ M each
Template DNA	x μ g	x μ g	/
HS Taq DNA Polymerase (5 U/ μ L)	0.4 μ L	1 μ L	0.1 U/ μ L
Nuclease-free Water	To 20 μ L	To 50 μ L	/

Note: Mix thoroughly before use.

- 1) Recommended template amounts (50 μ L reaction): Mammalian genomic DNA: 0.1 ~ 1 μ g; *E. coli* genomic DNA: 10 ~ 100 ng; Plasmid DNA: 0.1 ~ 10 ng
- 2) Enzyme volume can be adjusted from **0.25 ~ 1 μ L**. Higher volume usually improves yield but may reduce specificity.
- 3) Optimal final Mg²⁺ concentration: **1.5 ~ 2 mM**. Titrate with 25 mM MgCl₂ at 0.2 ~ 0.5 mM increments if needed.

2. PCR Protocol

Step	Temp		Time	Cycles
	eratu	re		
Initial denaturation	95°C		5 min	1
Denaturation	95°C		30 sec	30~35
Annealing	55°C		30 sec	
Extension	72°C		30 s / kb	
Final extension	72°C		10 min	1

Notes:

- 1) Extend initial denaturation to **10 min** for complex templates.
- 2) Annealing temperature: adjust based on primer T_m; typically **2 ~ 5 °C below T_m**.
- 3) Extension time: 30 sec / kb recommended; **60 sec / kb** for gene cloning to maximize yield.